



Edwards Lifesciences




Edwards Lifesciences



Hemodynamic Monitoring

Components - Heart



Hemodynamic Monitoring.

- **1) Pressure**

Arterial pressure , Central Venous Pressure

Pulmonary Artery Pressure (PAP) ,

Pulmonary Artery Wedge Pressure (PAWP)

Pressure

- **Catheter**
- **Pressure Transducer**
- **Monitor / Pressure module.**

Catheters

- **Arterial catheter (AP)**
- **Central Venous Catheter (CVC).....CVP**
- **Pulmonary Artery Catheter (PA Catheter)**
(CVP , PAP, PAWP)
- **Thermodilution Catheter. (CO Catheter)**
(CVP, PAP,PAWP,CO)

An anatomical diagram of the heart and pulmonary circulation. The heart is shown in a cross-section, with the right and left ventricles and atria. The pulmonary artery is shown branching out to the right lung, and the pulmonary veins are shown returning blood to the left atrium. The diagram is overlaid with a list of hemodynamic parameters in bold, italicized white text.

RIGHT ATRIAL PRESSURE
RIGHT VENTRICULAR PRESSURE
PULMONARY ARTERY PRESSURE
PULMONARY ARTERY WEDGE PRESSURE
CARDIAC OUTPUT
SYSTEMIC VASCULAR RESISTANCE
PULMONARY VASCULAR RESISTANCE
STROKE INDEX
VENTRICULAR STROKE WORK INDEX
CORONARY PERFUSION PRESSURE

Multi-Med

Central Venous Catheters (CVC) - CVP



- **16 cm & 20 cm lengths**
- **Single lumen to quad lumen (1L – 4L)**
- **Outer diameter (3 French – 8.5 French)**
- **Heparin coating (AMC THROMBOSHIELD)**



Edwards Lifesciences

Multimed CVC

Double Lumens

Distal /Brown 16 G..... CentralCVP

Proximal/white 16 G ...5cm...Medicine /TPN/Blood Sampling

Triple Lumens

.Distal/Brown 16G... Central...CVP measurement

Proximal /White 18G... 5 cm ..Blood Sampling/Medicine

Medial/Blue 18G.... 2 cm...Medicine

Multimed Single Lumen

- | Catheter Diameter | Volume | Flow Rate |
|----------------------------|----------------|-------------------|
| Distal Tip 14 G /6F | 0.40 CC | 6016 ml/hr |
| Distal Tip 16 G/5F | 0.23 CC | 3280 ml/hr |
| Distal Tip 20 G/3F | 0.10 CC | 1260 ml/hr |
- **Length 20 cm**

Multimed Double Lumens

- **Double Lumen (Flow rate)**

16 cm

20 cm

- **Proximal /16G 3,620 ml/hr 3,292 ml/hr**

(Exit port 5 cm)

- **Distal /16G 3,608 ml/hr 3,200 ml/hr**

(Exit port at the Tip)

Multimed Triple Lumens

- **Triple Lumens**

16 cm

20 cm

- **Proximal/18G**

1,670 ml/hr

1,420 ml/hr

(Exit port at 5 cm)

- **Medial/18G**

1,500 ml/hr

1,300 ml/hr

(Exit port at 2 cm)

- **Distal/16G**

3,510 ml/hr

3,160 ml/hr

(Exit port at the Tip)

Central Venous Catheter

- **1) High Flow Central Venous Catheter**

Infuse Fluids up to
50% Faster with *High-Flow*
Central Venous Catheters



CardioVascular Group
Edwards Critical-Care Division



Edwards Lifesciences

Central Venous Catheter

- **2) Coated Central Venous Catheter**
(Antimicrobial + Heparin) , Double protection
Benzalkonium Chloride + Heparin


**Some brand only - Sulfadiazine , Silver Nitrate ,
No Heparin , Only one side.**

Central Venous Catheter

3) Central Venous Catheter with Introducer capability

AVA

Advanced Venous Access Devices




AVA HF

Advance Venous Access

AVA Devices

Competencies in extrusion
Innovation and technology

Flexible Inner Lumens



Increased flow rates

AVA HF (High-Flow)

Advanced Venous Access Devices

Make sure that AVA HF is covering your

H♥RTTS

High **R**isk ♥ (cardiac) **T**rauma

Re-do CABG or Valve
CABG (multi-vessel)
Unstable patient

Accident

GSW

All Trauma centers

Transplant

Liver
Heart
Lung

Surgical

Abdominal Aortic

Aneurysm

Vascular

Anticipate difficulty (prostate)




Edwards Lifesciences

AVA 3Xi

Advanced Venous Access Devices

Central Venous Catheter with Introducer capability



Vantex
WITH OLIGON


**Revolutionizing the Science
of Antimicrobial Protection**




First CVC with an integrated antimicrobial material


OLIGON

- Antimicrobial compound
 - *Silver*, platinum, and carbon
- Body fluids interact with silver and platinum particles in the material, causing a release of silver ions





PreSep[®]
Central Venous Oximetry
ScvO₂



Edwards PreSep Central Venous ScvO₂ Oximetry Catheter

- Up to 50% of patients resuscitated from shock may have continued global tissue hypoxia (i.e., increased lactate and decreased ScvO₂) even with the normalization of vital signs and central venous pressure¹
- Reduced central blood volume is reflected more clearly with ScvO₂ than in CVP²
- ScvO₂ saturation is a reliable and sensitive method for detecting blood loss³




PreSep Central Venous Oximetry Catheter* Specifications:

| Model Number | Lumens | Length (cm) | Size F (mm) | Continuous ScvO ₂ | Lumen Size Gauge (mm) | | | Recommended Dilator F (mm) | Minimum Guidewire Size Inch (mm) | AMC Thrombosshield** |
|--------------|--------|-------------|-------------|------------------------------|-----------------------|-----------|-----------|----------------------------|----------------------------------|----------------------|
| | | | | | Distal | Proximal | Medial | | | |
| X3820HK | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |
| X3820K | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |
| X3820HS*** | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |

*PreSep catheters are designed for use with Edwards Lifesciences SAT-2 device, Explorer monitor, Vigilance monitor and OM2 optics module to continuously monitor ScvO₂.

**All model numbers with an "H" contain AMC Thrombosshield, an antibacterial heparin coating which decreases viable microbe count on surface of product during handling and placement.

*** Model numbers with an "S" do not contain xylocaine local anesthesia.



ScvO₂ a sensitive indicator of changes in:

Oxygenation:

FiO₂ Ventilation

Cardiac Output:


Heart Rate, Preload, Afterload, Contractility

Hemoglobin

Bleeding, Hemodilution

Metabolic Demand

Shivering, Work of breathing, Fever, Seizures



Compatible with Edwards Swan-Ganz SvO₂ Optics modules and computers

- Vigilance monitor
- Explorer monitor
- SAT-2 device



Soft Tip. Helps reduce the likelihood of complications resulting from vessel perforation.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards Lifesciences, Edwards, the stylized E logo, PreSep and SAT-2 are trademarks of Edwards Lifesciences Corporation. AMC Thrombosshield, Explorer, Swan-Ganz and Vigilance are trademarks of Edwards Lifesciences Corporation and are registered in the U.S. Patent and Trademark office.

© 2004 Edwards Lifesciences LLC
All rights reserved. ARXXXX


References

1. Rivers M, et al. Central venous oxygen saturation monitoring in the critically ill patient. *Curr Opin Crit Care* 2001; 7(3):204-11
2. Madsen P et al., Central venous oxygen saturation during hypovolaemic shock in humans. *Scand J Clin Lab Invest* 1993; 53:67-72
3. Scalea TM et al., Central venous oxygen saturation: a useful clinical tool in trauma patients. *J Trauma* 1990; 30:1539-1543




Edwards Lifesciences


ScvO₂ measurement



The Swan-Ganz®

Catheter





Right Atrium




“x” = Atrial Diastole


“y” = Atrial Emptying

Right Atrial Pressure

Mean 2 - 6 mm Hg



Right Ventricle





- 1. Isovolumic Contraction**
- 2. Rapid Ejection**
- 3. Reduced Ejection**
- 4. Isovolumic Relaxation**
- 5. Early Diastole**
- 6. Atrial Systole**
- 7. End-Diastole**


Right Ventricular Systolic Pressure


15 - 25 mm Hg

Right Ventricular End-Diastolic Pressure

0 - 8 mm Hg








Pulmonary Artery Systolic Pressure


15 - 25 mm Hg

Pulmonary Artery Diastolic Pressure


8 - 15 mm Hg

Pulmonary Artery





Pulmonary Artery Wedge



“a” Wave = Atrial Contraction

“x” Descent = Atrial Diastole


“v” Wave = Passive Atrial Filling



“y” Descent = Atrial Emptying

Pulmonary Artery Wedge Pressure

Mean

6 - 12 mm Hg







The image shows a TruWave Disposable Pressure Transducer, which is a clear plastic device with a blue cable attached. The cable has a blue connector with a Snap-Tab. The device is shown in a close-up view, highlighting its ergonomic design and the flow-through design.

TruWave™ Disposable Pressure Transducer

An Easy, Convenient System...

The ease of using the TruWave™ Disposable Pressure Transducer is immediately apparent – from the ergonomical design that fits your hand, to the cable connector* that has been designed to connect and disconnect easier. The unique flow-through design provides an unimpeded fluid path for easy filling of the system. The protective connector sheath is water resistant and connects without pins. The test port minimizes trouble shooting by testing both the monitoring cables and disposable pressure transducer.

An Accurate System...

The backside test port lets you verify the accuracy of the system quickly and easily. The convenient Snap-Tab™ flush device, designed for a sure grip, can be pulled from any

direction to flush your system and will generate a square-wave test pattern. With greater Snap-Tab™ sensitivity, you “feel” how fast you’re flushing the system.

And With VAMP® Kits, A Complete System. Paired with VAMP® Kits (Venous Arterial blood Management Protection system), you have a complete and reliable closed, needleless blood sampling system with accurate pressure readings in one complementary unit.

For years, the Edwards Critical-Care Division of Baxter Healthcare has brought you the benefit of its critical care experience. You told us what you wanted, and we listened.

Edwards Critical-Care Division

Baxter

*Patent Pending



Edwards Lifesciences

Hemodynamic

- **2) Volume / Cardiac output**
Swan Ganz Thermodilution
- **3) Oxgenation**
Swan Ganz with SVO2


Hemodynamic Monitoring

- **2) Volume**
 - **Stroke Volume (SV)**
 - **Ejection Fraction (EF)**
 - **Continuous end diastolic volume (CEDV)**
 - **End Systolic Volume (ESV)**
 - **Cardiac output (CO, CCO)**

Hemodynamic Monitoring

- **3) Oxygen Profile**
 - **Delivery (DO₂) , Consumption (VO₂)**
 - **Mixed Venous Oxygen Saturation (SVO₂)**
 - **Central Venous Oxygen Saturation (ScVO₂)**

Hemodynamic Measurements



Cardiac Output (CO)

$$\text{CO} = \text{Heart Rate (HR)} \times \text{Stroke Volume (SV)}$$



Cardiac Index (CI)

$$CI = \frac{CO}{\text{Body Surface Area (BSA)}}$$



Bois, E.F. Basal Metabolism
Health Disease, Lea and
Finger, 1936.




CO = 4 - 8 liters/min


CI = 2.5 - 4.0 liters/min/m²



Thermodilution

General Principles





Technical Factors Necessary for Accurate Bolus CO

- **Appropriate technique**
- **Appropriate catheter position**
- **Accurate temperature of the injectate**
- **Accurate volume of injectate**
- **Appropriate computation constant**
- **Consistent timing of injection**
- **Consistent averaging strategy**

Appropriate Technique

- **Steady, consistent injections**
- **10 cc / 4 seconds**

Appropriate Catheter Position

- **Appropriate wedge tracing**
- **Balloon inflation volume: 1.25 - 1.5 cc**
- **Appropriate RA tracing**

Accurate Temperature of the Injectate

- **1°C error in injectate temperature:**

2.7 % using 0°C (iced) injectate

**7.7% using 24°C (room temperature)
injectate**

Levett J & Replogle RL; 1979

Accurate Injectate Volume

- **A 0.5 ml error in a 5 ml injection will cause a 10% error.**

Weissman, (1987) Measuring Oxygen Uptake.

Consistent Timing of Injection

- **Issue of accuracy vs reproducibility**
- **Reproducibility: variation from respiratory cycle reduced.**
- **If CO determinations shot throughout the cycle: variance reported to be as high as 70%.**

Jansen, 1981.

Case Studies: What has changed in this patient?

| | | |
|--------------|---------------|---------------|
| • C.O./C.I | 2.4 / 1.50 | 4.91 / 3.07 |
| • PAWP | 13 | 14 |
| • PAP (S/D) | 36 / 18 | 40 / 15 |
| • BP (S/D/M) | 106 / 68 / 81 | 110 / 70 / 85 |
| • HR | 100 | 90 |
| • RAP | 6 | 7 |
| • SVR | 2643 | 1270 |

Computation Constant: Function of lumen size, injectate volume, injectate temperature


- **During the night shift, the RN increased the volume of injectate from 5 to 10 cc to enhance reproducibility.**
- **CO Wrong x CC Right / CC Wrong**
- **$2.4 \times 0.561 / 0.259 = 5.2$**

Consistent Averaging Strategy


- **Delete values associated with poor curves, alerts.**
- **Attain triplicate measurements.**
- **Common practice: delete values outside a median value of 10%.**

CONTINUOUS CARDIAC OUTPUT


Continuous Cardiac Output –Vigilance I



Vigilance II




Swan Ganz –CCO Catheter



- **Assess
Right Atrial Tracing**


- **Assess PAWP**
- **Assess Balloon
Inflation Volume**



Input Signal

- **Small energy signals (indicator) are infused directly into the blood in an apparently random (but actually repeating), on-off pattern.**


Random, on-off pattern




Output Signal

- **The resulting blood temperature changes are detected at the thermistor in the pulmonary artery.**



Pulmonary artery temperature changes



Continuous Cardiac Output Washout Curve



Cross Correlator



Stroke Volume

Preload

Afterload


Contractility



Swan Ganz –CCO /Svo2 /CEDV /EF

- **Continuous cardiac output**
- **Mixed Venous oxygen saturation**
- **Continuous End diastolic volume**
- **Ejection Fraction**

Swan Ganz – CCO/SVO2





FloTrac Sensor
&
Vigileo Monitor

*A New Platform for Minimally Invasive
Hemodynamic Monitoring*

Preecha Bhandtvej

Minimally Invasive Monitoring


Vigileo Monitor



FloTrac Sensor





FloTrac: Easy, Accurate and Convenient



- **From a standard arterial line**
- **Monitors an important subset of parameters**
 - Cardiac Output, Oxygenation
- **Applicable to many more patients**


System Configuration

PreSEP Catheter
(central vein)




Oximetry

Vigileo
Instrument



**Cardiac
Output**

FloTrac Sensor
(peripheral artery)



Edwards PreSep Central Venous ScvO₂ Oximetry Catheter

- Up to 50% of patients resuscitated from shock may have continued global tissue hypoxia (i.e., increased lactate and decreased ScvO₂) even with the normalization of vital signs and central venous pressure¹
- Reduced central blood volume is reflected more clearly with ScvO₂ than in CVP²
- ScvO₂ saturation is a reliable and sensitive method for detecting blood loss³




PreSep Central Venous Oximetry Catheter* Specifications:

| Model Number | Lumens | Length (cm) | Size F (mm) | Continuous ScvO ₂ | Lumen Size Gauge (mm) | | | Recommended Dilator F (mm) | Minimum Guidewire Size Inch (mm) | AMC Thrombosshield** |
|--------------|--------|-------------|-------------|------------------------------|-----------------------|-----------|-----------|----------------------------|----------------------------------|----------------------|
| | | | | | Distal | Proximal | Medial | | | |
| X3820HK | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |
| X3820K | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |
| X3820HS*** | 3 | 20 | 8.5 (2.83) | • | 15 (1.77) | 18 (1.33) | 18 (1.33) | 10.5 (3.5) | 0.32 (0.8) | • |

*PreSep catheters are designed for use with Edwards Lifesciences SAT-2 device, Explorer monitor, Vigilance monitor and OM2 optics module to continuously monitor ScvO₂.

**All model numbers with an "H" contain AMC Thrombosshield, an antibacterial heparin coating which decreases viable microbe count on surface of product during handling and placement.

*** Model numbers with an "S" do not contain xylocaine local anesthesia.



ScvO₂ a sensitive indicator of changes in:

Oxygenation:

FiO₂ Ventilation

Cardiac Output:


Heart Rate, Preload, Afterload, Contractility

Hemoglobin

Bleeding, Hemodilution

Metabolic Demand

Shivering, Work of breathing, Fever, Seizures



Compatible with Edwards Swan-Ganz SvO₂ Optics modules and computers

- Vigilance monitor
- Explorer monitor
- SAT-2 device



Soft Tip. Helps reduce the likelihood of complications resulting from vessel perforation.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards Lifesciences, Edwards, the stylized E logo, PreSep and SAT-2 are trademarks of Edwards Lifesciences Corporation. AMC Thrombosshield, Explorer, Swan-Ganz and Vigilance are trademarks of Edwards Lifesciences Corporation and are registered in the U.S. Patent and Trademark office.

© 2004 Edwards Lifesciences LLC
All rights reserved. ARXXXX


References

1. Rivers M, et al. Central venous oxygen saturation monitoring in the critically ill patient. *Curr Opin Crit Care* 2001; 7(3):204-11
2. Madsen P et al., Central venous oxygen saturation during hypovolaemic shock in humans. *Scand J Clin Lab Invest* 1993; 53:67-72
3. Scalea TM et al., Central venous oxygen saturation: a useful clinical tool in trauma patients. *J Trauma* 1990; 30:1539-1543



Edwards Lifesciences

ScvO₂ measurement



Sepsis, EGDT

Sepsis: Find it Fast.
Manage it Early.




Edwards PreSep™ Central Venous Oximetry Catheter



Early Goal-Directed Therapy (EGDT)
with PreSep yields significant
reductions in sepsis related mortality.

We observed a
lower mortality
rate in patients
with septic
shock assigned
to early
goal-directed
therapy
(42.5 percent)
than in those
assigned to
standard therapy
(56.8 percent)²



Results:²


34% Reduction of in-hospital mortality

3.8 Reduction in hospital days

\$12,000 Reduction in total hospital charges

Sepsis, EGDT

Case Study



History of Present Illness:

85-year-old male with a history of hypertension, atherosclerotic heart disease, and congestive heart failure. He presents with cough and shortness of breath that started while making breakfast. He became ill while visiting his wife in the hospital and decided to come to the Emergency Department from the hospital floors.

Physical Examination:

Temperature 34.5, pulse 84, blood pressure 150/98, respiratory rate 28. Elderly male who was awake, talkative, alert and in no acute distress. HEENT: mild jugular venous distention at 30 degrees, dry mucous membranes. LUNGS: crackles in the left lung base. CARDIOVASCULAR: Regular rate and rhythm. No auscultatory murmur. ABDOMEN: Mildly distended and tympanitic below the umbilicus. RECTAL: Heme-negative, brown stool with very enlarged prostate.

Diagnosis:

1. Urosepsis secondary proteus mirabilis and obstructive uropathy.
2. Pneumonia – pseudomonas aeruginosa.
3. Mild Consumptive coagulopathy.
4. Decompensated heart failure-exacerbated by sepsis.
5. Acute renal insufficiency secondary to obstructive uropathy and acute tubular necrosis.

| Hour | HR | MAP | CVP | RR | Urine Output | HGB | ScvO ₂ | SaO ₂ | Lactic Acid | Diagnosis and Medical Intervention |
|------|----|-----|-----|----|--------------|------|-------------------|------------------|-------------|--|
| 0 | 79 | 110 | 33 | 32 | 0 | 11.5 | 98.7 | 14 | 5.3 | Chest x-ray positive for left lower lobe pneumonia. Foley catheter site erythematous but unremarkable, urinalysis consultation, increased CVP and pro-renal azotemia indicates myocardial insufficiency. |
| 1 | 88 | 83 | 28 | 30 | 500 | | 99 | 28 | | Suprapubic catheterization performed, large amount of cloudy urine is present (post-obstructive diuresis), urinalysis positive for urinary tract infection. Antibiotics started for pneumonia and urinary tract infection. Dobutamine started at 2.5 mcg/kg/min* |
| 2 | 84 | 90 | 26 | 32 | 30 | | 97 | 49 | | Dobutamine increased to 5 mcg/kg/min |
| 3 | 83 | 89 | 20 | 28 | 25 | | 98 | 29 | 4.1 | Dobutamine increased to 7.5 mcg/kg/min. Sudden drop in ScvO ₂ , patient became agitated, tachycardic which necessitated intubation and mechanical ventilation. |
| 4 | 88 | 87 | 15 | 30 | 20 | | 98 | 64 | | Patient sedated and dobutamine increased to 10 mcg/kg/min |
| 5 | 90 | 79 | 6 | 16 | 35 | | 97 | 69 | 3.6 | Fluid challenge given and dobutamine increased to 12.5 mcg/kg/min |
| 6 | 97 | 79 | 10 | 16 | 55 | | 99 | 75 | 2.9 | Repeat fluid challenge and additional sedation provided |
| 7 | 95 | 85 | 15 | 16 | 75 | 10 | 99 | 82 | 1.5 | Patient transported to the intensive care unit, sedated 3 days later and discharged 7 days later to home |

Sepsis patients may be at high risk of irreversible organ failure, even in the presence of normal vital signs.

"Up to 50% of patients resuscitated from shock may have continued global tissue hypoxia (i.e. increased lactate and decreased ScvO₂) even with the normalization of vital signs and central venous pressure."³

Protocol for Early Goal-Directed Therapy for Sepsis.

"The early identification of patients with insidious illness (global tissue hypoxia accompanied by stable vital signs) makes possible the early implementation of goal directed therapy...we conclude that goal-directed therapy provided at the earliest stages of severe sepsis and septic shock... has significant short-term and long-term benefits."¹

Screen Early for At-Risk Patients:


2 Signs of the Systemic Inflammatory Response Syndrome (SIRS) +


Sign of Global Tissue Hypoxia

Temp C° < 36°C or > 38°
HR > 90 beats per/m
Resp > 20 breaths per/m
WBC > 12.0 x 10⁹/liter (L)

SBP < 90 mm Hg
or
Lactate 4 mmol/L

Early Treatment:





Parameters

| Label | Parameter | Range/Units |
|----------------------|------------------------------------|--|
| CO | Arterial Pressure Cardiac Output | 1.0 - 20.0 L/min |
| ScvO ₂ ** | Central Venous Oxygen Saturation | 0 - 99% |
| SvO ₂ ** | Mixed Venous Oxygen Saturation | 0 - 99% |
| CI | Cardiac Index | 0 - 20.0 L/min/m ² |
| SV | Stroke Volume | 0 - 300 ml/beat |
| SVI | Stroke Volume Index | 0 - 200 ml/beat m ² |
| SVV | Stroke Volume Variation | 0 - 99% |
| SVR | Systemic Vascular Resistance | 0 - 3,000 dynes-sec/cm ⁵ (0 - 300.0 kPa-sec/l) |
| SVRI | Systemic Vascular Resistance Index | 0 - 6,000 dynes-sec-m ² /cm ⁵ (0 - 600.0 kPa-sec-m ² /l) |





Edwards

Helping patients is our life's work, and

life is now



Edwards Lifesciences